

<i>ECG Functions</i>	<b>Green v.05.101</b>	<b>Green v.06.101</b>	<b>Orange</b>	<b>Grey v.07.205</b>	<b>Grey v.07.305</b>	<b>Gold3 v.202</b>	<b>Gold3 v.305</b>	<b>Coral</b>	<b>Amber</b>
Display	graphic 4'	colour LCD TFT 4,3" touch screen	graphic 4'	colour LCD TFT 7" touch screen	colour LCD TFT 7" touch screen	colour LCD TFT 10,4" touch screen	colour LCD TFT 10,4" touch screen	-	-
Keyboard	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	membrane alphanumeric with functional keys	-	-
Communication interface	micro USB	micro USB	-	Ethernet, USB device, USB host	Wi-Fi, Ethernet, USB device, USB host	Ethernet, USB device, USB host, PS/2	Wi-Fi, Ethernet, USB device, USB host, PS/2	USB 1.1	wireless transmitter AsLINK
Dimensions (mm)	220x153x55	220x153x55	258x199x50	258x199x50	258x199x50	370x372x94	370x372x94	118x79x27	130x96x30
Weight (kg)	<0,6	<0,6	<1,3	<1,3	<1,3	<4,3	<4,3	<0,2	<0,25
<b>Printout</b>									
Paper (roll / mm)	58	58	112	112	112	210	210	-	-
External printer	yes (from PC)	yes (from PC)	no	yes (from ECG or PC)	yes (from ECG or PC)	yes (from ECG or PC)	yes (from ECG or PC)	yes (from PC)	yes (from PC)
<b>ECG signal</b>									
Displayed leads	1/3/6/12	1/3/6/12	3/6/12	1/3/6/12	1/3/6/12	1/3/6/12	1/3/6/12	3/6/12	3/6/12
Printed leads	1/3/6/12	1/3/6/12	3/6/12	1/3/6/12	1/3/6/12	1/3/6/12	1/3/6/12	3/6/12	3/6/12
Working modes	auto, manual	auto, manual	auto, manual	auto, manual, automanual	auto, manual, automanual	auto, manual	auto, manual, automanual	auto, manual	auto, manual
ECG back recording	no	no	no	yes	yes	no	yes	no	no
Rhythms' printout	no	no	no	yes	yes	no	yes	no	no
LONG mode	no	no	no	yes	yes	no	yes	no	no
Sampling frequency (Hz)	8000 Hz	8000 Hz	2000 Hz	8000 Hz	8000 Hz	8000 Hz	8000 Hz	1600 Hz	500 Hz
Sensitivity (mm/mV)	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20	2,5/5/10/20
Speed (mm/s)	5/10/25/50	5/10/25/50	5/10/25/50	5/6,25/10/12,5/25/50	5/6,25/10/12,5/25/50	5/ 6,25/ 10/ 12,5/ 25/ 50/ 100	5/ 6,25/ 10/ 12,5/ 25/ 50/ 100	5/10/25/50	5/10/25/50
Internal memory (examinations)	5-1000	5-1000	last automatic examination	1000	1000	2000	2000	PC drive	PC drive
Work in a battery mode	up to 130 automatic examinations	up to 130 automatic examinations	up to 130 automatic examinations	up to 130 automatic examinations	up to 130 automatic examinations	up to 300 automatic examinations	up to 300 automatic examinations	no	up to 12 hours
<b>Other</b>									

Arkusz1

Analysis and interpretation	yes	yes	yes	yes	yes	yes	yes	yes	yes
Spirometry screening	no	no	no	yes	yes	yes	yes	no	no
LAN / Internet	no	no	no	yes (cable)	yes (Wi-Fi, cable)	yes (cable)	yes (Wi-Fi, cable)	no	no
HL7	no	no	no	no	yes	no	yes	no	no
ECG-Mail	no	no	no	yes	yes	yes	yes	yes	yes
Electronic Medical Records (EMR)	no	no	no	yes	yes	yes	yes	yes	yes
INOP	yes	yes	yes	yes	yes	yes	yes	yes	yes